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Irrigation solution, sterile and endotoxin-free

1,000 mL of solution contains:	
Sodium chloride	6.00 g
Potassium chloride	0.30 g
Calcium chloride 2H <sub>2</sub> O	0.20 g
Magnesium chloride 6H <sub>2</sub> O	0.20 g
Sodium lactate	3.1 g
in water for injections	
<b>Electrolytes:</b>	<b>mmol/L</b>
Sodium	130
Potassium	4.02
Calcium	1.36
Magnesium	0.98
Chloride	111
Lactate	27.7

**Therapeutic indications**

For rinsing and cleansing during surgical intervention  
For irrigation of wounds and burns  
For irrigation and cleansing in stoma care  
For irrigation during endoscopic investigations in body cavities  
For intraoperative irrigation during arthroscopy with mechanical instruments  
For postoperative irrigation  
For mechanical irrigation of the eyes (the use of 30 mL container is recommended)

**Posology**  
The amount of fluid depends on the clinical condition.

**Overdose**  
Overdosage during wound irrigation and intraoperative application can lead to absorption of irrigation fluid, that can lead to hypervolaemia in the form of hypotonic hyperhydration (hyponatraemia). Headache, nausea, restlessness and disorientation are symptoms of hyperhydration. In severe cases, coma-like states can occur.

**Countermeasures in case of overdose**  
Discontinuation of irrigation, correction of the water and electrolyte balance depending on the particular clinical situation.

**Warnings and precautions**  
Do not use for infusion. Only use the solution if the seal is undamaged and the solution is clear. Opened containers must not be saved, quantities left over must be discarded. Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death

## Ringer Lactate B. Braun

of the patient. Keep out of the sight and reach of children. If there is a risk that large quantities of solution can be ingested, as is the case when using in large body cavities or on large wound areas, the solution must not be used on patients in a state of hyperhydration. Caution is advised in the following situations:

– Hypertonic dehydration,  
– Hyperkalaemia,  
– Hypermtraemia,  
– Hyperchloraemia,  
– Renal insufficiency with a tendency to hyperkalaemia,  
– Disorders in which reduced sodium intake is indicated e.g. heart failure, generalised oedema, pulmonary oedema, hypertension, eclampsia, severe renal insufficiency.

Monitoring should include checks for water and serum electrolyte balance. Ringer Lactate B. Braun must not be used for intrathecal irrigation procedures.  
**Undesirable effects**  
None are known.  
**Expiry date**  
The solution must not be used after the expiry date which is stated on the label.  
**Date of last revision**  
0310

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**B. Braun Melsungen AG**