


# AESCULAP® Detachable Rongeurs



# Noir<sup>®</sup> Coated Detachable Rongeurs

## Cleanability

Ever increasing demands in the reprocessing of surgical instruments require professional and efficient processes, together with intelligent product solutions, which meet modern hygiene standards.




## Noir<sup>®</sup> coating

All models come with the proven Noir<sup>®</sup> coating for the reduction of distracting light reflections.<sup>1</sup>

## The choice is yours

86 models with working lengths from 150 mm to 280 mm, jaw widths from 1.5 mm to 5 mm, standard and semi-serrated jaws in straight and 150° upward curved.



<sup>1</sup> The gloss level of a Noir<sup>®</sup>-coated instrument was compared to the gloss level of an uncoated instrument by means of a glossmeter. The reflections (gloss levels were determined at 6 different measure points) were reduced by 60% (crosswise) and 80% (lengthwise) on average.



## Handling

The golden button indicates the detachability. In the open cleaning position all parts remain connected.


## Identification

- Working length
- Jaw width

## Disassembly steps

- 1 Push the golden button
- 2 Pull down the front handle
- 3 Lift the slider

# Noir<sup>®</sup> Coated Detachable Rongeurs



## Rongeurs standard



Working length	Jaw	1.5 mm	2 mm	3 mm	4 mm	5 mm
150 mm		FH801B	FH802B	FH803B	FH804B	FH805B
150 mm		FH901B	FH902B	FH903B	FH904B	FH905B
180 mm		FH811B	FH812B	FH813B	FH814B	FH815B
180 mm		FH911B	FH912B	FH913B	FH914B	FH915B
200 mm		FH821B	FH822B	FH823B	FH824B	FH825B
200 mm		FH921B	FH922B	FH923B	FH924B	FH925B
230 mm		FH891B	FH892B	FH893B	FH894B	FH895B
230 mm		FH931B	FH932B	FH933B	FH934B	FH935B
280 mm		-	-	FH837B	FH838B	FH839B




### Rongeurs semi-serrated



Working length	Jaw	1.5 mm	2 mm	3 mm	4 mm	5 mm
150 mm		FH851B	FH852B	FH853B	FH854B	FH855B
150 mm		FH951B	FH952B	FH953B	FH954B	FH955B
180 mm		FH861B	FH862B	FH863B	FH864B	FH865B
180 mm		FH961B	FH962B	FH963B	FH964B	FH965B
200 mm		FH871B	FH872B	FH873B	FH874B	FH875B
200 mm		FH971B	FH972B	FH973B	FH974B	FH975B
230 mm		FH881B	FH882B	FH883B	FH884B	FH885B
230 mm		FH981B	FH982B	FH983B	FH984B	FH985B
280 mm		-	-	FH847B	FH848B	FH849B

# Instrument Organisation System






## AESCULAP® Instrument Organisation System

Surgical instruments are valuable assets for every hospital, which is why they should be processed as efficiently and gently as possible.

Therefore, protecting the instruments during sterilisation, as well as in transport to sterile supply stores and to the operating room, is of crucial importance.

The new Aesculap organisation system for instrument trays facilitates the storage especially of delicate instruments and small items.



The main product trademark "AESCULAP" and the product trademark "Noir" are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

## AESCULAP® – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | [www.bbraun.com](http://www.bbraun.com)